



# LAGM Laser Ablation GED\_MSAG System




IAS Inc.

A 300 mm wafer sample can be analyzed by Laser Ablation (LA) ICP-MS without using a small, enclosed cell. Generated particles by LA are aspirated by an ejector and introduced to ICP-MS via GED.

## Features


- ◆ Femto-second laser and Galvanometer mirror ablates up to 12" (300 mm) wafer precisely.
- ◆ A wafer sample is laser ablated without using a small, enclosed cell, and the generated particles are aspirated by the ejector and introduced to the ICP-MS via Gas Exchange Device (GED).
- ◆ Dual syringe model MSAG\_DS enables quantitative analysis under a dry plasma condition of ICP-MS using the method of standard addition.
- ◆ Ultra-trace level of spot and wafer bevel analysis is available.
- ◆ Deeper level, > 300  $\mu\text{m}$ , of profile and impurity analysis is available.
- ◆ Two load ports from 12" (300 mm) to 6" (150 mm) wafer, aligner and wafer transfer robot are integrated for fully automated operation (Option)
- ◆ Alternative sampling stage is available for non-wafer samplers.



Schematic diagram of LA-GED-MSAG-ICP-MS

# Analysis of spiked Si wafer

- TRA analysis of spiked Si wafer sample
- 5 uL of 0.1 ppb mixed standard was spiked on a Si wafer and dried.



LA frequency: 10 kHz, Galvanometer mirror scan speed: 50 mm/sec, Galvanometer mirror scan width: 5 mm (Y axis direction),  
 Wafer stage movement speed: 0.0476 mm/sec (X axis), Wafer stage movement distance : 40 mm (X axis)  
 ICP-MS data acquisition mode: Transient analysis, Dwell time: 50 msec/mass, Date acquisition interval : 1.62 sec


- Quantitative analysis of SiC wafer using the method of standard addition

Analyte	Mass	Ar						SiC						Subtraction (ag)	Conc. (wt. ppb)
		MSAG 10 ppb STD flow (uL/min)			Corr.	ag/count	ag	MSAG 10 ppb STD flow (uL/min)			Corr.	ag/count	ag		
		0	1.5	3				0	1.5	3					
Li	7	16	174,989	367,620	1.000	1.4	22	12	176,033	367,165	1.000	1.4	16	(6)	< 51
Na	23	267	259,241	501,581	1.000	1.0	267	154	259,383	496,774	1.000	1.0	155	(111)	< 137
Mg	24	3	36,759	76,255	1.000	6.6	18	30	42,229	77,791	0.999	6.4	195	177	647
Al	27	469	196,939	400,024	1.000	1.3	587	4,853	196,846	398,317	1.000	1.3	6167	5581	20,400
K	39	855	89,725	193,945	0.999	2.6	2214	848	88,370	196,667	0.998	2.6	2165	(49)	< 814
Ca	40	12	77,147	165,683	0.999	3.0	37	6,095	78,132	171,641	0.997	3.0	18409	18372	67,000
Ti	48	4	3,684	7,533	1.000	66.4	266	30	4,401	7,658	0.996	65.5	1988	1723	6,280
V	51	7	42,665	93,469	0.999	5.3	37	9	42,062	96,303	0.997	5.2	48	11	< 174
Cr	52	2,184	98,824	204,904	1.000	2.5	5386	2,301	97,829	209,541	0.999	2.4	5551	165	< 1270
Mn	55	241	195,973	406,693	1.000	1.2	297	149	193,386	417,034	0.999	1.2	179	(118)	< 160
Fe	56	508	158,139	333,800	1.000	1.5	762	4,205	210,726	347,230	0.993	1.5	6130	5368	19,600
Fe	54	113	12,092	25,652	0.999	19.6	2212	282	14,793	26,556	0.998	19.0	5373	3161	11,500
Co	59	28	112,862	243,015	0.999	2.1	58	24	111,429	248,375	0.998	2.0	48	(9)	< 108
Ni	60	77	27,434	58,610	0.999	8.5	655	93	27,154	60,136	0.998	8.3	772	116	< 877
Cu	63	649	78,263	165,413	0.999	3.0	1970	986	77,793	170,014	0.999	3.0	2917	946	3,450
Cu	65	304	35,324	75,103	0.999	6.7	2032	470	34,998	76,734	0.999	6.6	3084	1052	3,840
Zn	66	32	34,001	71,822	1.000	7.0	221	148	34,091	73,350	0.999	6.8	1009	788	2,870
Zn	64	12	54,735	115,774	1.000	4.3	52	191	54,951	119,118	0.999	4.2	804	753	2,750
Ga	71	508	127,800	266,371	1.000	1.9	955	555	130,726	271,534	1.000	1.8	1024	69	< 476
Ge	74	10	5,030	10,857	0.999	46.1	475	12	5,184	10,817	1.000	46.3	540	65	< 1,730
As	75	72	18,006	36,732	1.000	13.6	978	70	17,893	36,577	1.000	13.7	959	(19)	< 1,250
Rb	85	94	190,944	400,964	1.000	1.2	118	78	192,041	399,331	1.000	1.3	97	(20)	< 121
Sr	88	14	205,366	428,998	1.000	1.2	16	117	206,084	427,355	1.000	1.2	137	121	441
Zr	90	22	43,861	91,653	1.000	5.5	120	14	43,724	91,067	1.000	5.5	75	(45)	< 222
Mo	98	1,066	36,940	76,816	1.000	6.6	7034	1,116	36,398	78,275	0.999	6.5	7230	196	< 2,370
Ag	107	9	74,376	154,641	1.000	3.2	28	16	73,942	153,750	1.000	3.3	51	23	< 141
Cd	111	7	16,437	34,420	1.000	14.5	97	9	16,580	34,453	1.000	14.5	131	33	< 477
Sn	118	78	39,748	82,164	1.000	6.1	477	160	39,969	81,987	1.000	6.1	980	503	1,830
Sb	121	9	38,013	79,010	1.000	6.3	59	25	38,016	78,353	1.000	6.4	160	101	367
Cs	133	26	175,692	364,589	1.000	1.4	36	26	176,594	362,761	1.000	1.4	35	(0)	< 76
Ba	138	7	100,331	209,654	1.000	2.4	17	22	100,463	208,022	1.000	2.4	52	35	129
W	184	318	21,326	44,558	1.000	11.3	3591	323	21,279	43,969	1.000	11.5	3704	114	< 2,250
Pb	208	1	93,474	194,978	1.000	2.6	2	63	92,306	192,864	1.000	2.6	163	161	586
Pb	206	0	45,381	94,370	1.000	5.3	0	29	44,925	93,498	1.000	5.3	153	153	559
Si	28	9,346	15,582	21,355	1.000	4163.6	3.89.E+07	50,492	(*1) 57,201	(*2) 61429	0.992	4571.6	2.31.E+08	191,916,935	


LA frequency: 10 kHz, Galvanometer mirror scan speed: 50 mm/sec, Galvanometer mirror scan width: 1 mm (Y axis)  
 Wafer stage movement speed: 0.025 mm/sec (X axis),  
 ICP-MS data acquisition mode: Spectrum, Integration time: 1 sec/mass

# Analysis of 0.5 ppb 20 nm Au particles on Si wafer

- Single nano-particle analysis of Au particles on Si wafer sample.
- 3  $\mu\text{L}$  of 0.5 ppb Au particles standard solution in IPA was spiked on a Si wafer and dried.




X axis full scale : 5 sec



X axis full scale : 0.002 sec



LA frequency: 10 kHz, Galvanometer mirror scan speed: 50 mm/sec, Galvanometer mirror scan width: 15 mm (X Y axis direction), ICP-MS data acquisition mode: single nanoparticle mode, Dwell time: 0.1 msec/mass, Analysis time : 600 sec



Au sensitivity by MSAG : 2.94 ag/count  
 Average Au particle count : 28 counts  
 # of Au particle in 3  $\mu\text{L}$  : 19,083 pcs  
 # of Au particles detected : 16,009 pcs

Au nano-particle size calculated from MSAG was 20.1 nm  
 The number of particle recovery was 84%

# Analysis of GaN wafer




LA frequency: 10 kHz, Galvanometer mirror scan speed: 50 mm/sec, Galvanometer mirror scan width: 5 mm (Y axis direction), Wafer stage movement speed: 0.025 mm/sec (X axis), ICP-MS data acquisition mode: Single nano-particle analysis, Dwell time: 0.1 msec/mass, Analysis time : 60 sec




# Analysis of bevel


Photo of bevel area  
Laser ablation width : 250  $\mu\text{m}$




Analysis of Au particles on bevel  
(Two droplets of 20 nm Au particle solution)



Analysis of bevel of Si wafer (1st ablation): Several impurities were found at the same spot.



Analysis of spiked wafer (Two droplets of STD solution on the bevel area after 7 times cleaning):  
Two peaks were detected.



LAGM software communicates with ICP-MS software that enables fully automated operation.

- Calibration mode : Calibration curves while laser ablating a sample are obtained using *MSAG\_DS*.
- Full scan mode : Entire wafer is laser ablated according to the analysis time. LAGM software automatically calculates speed of X-Y-Z-θ stage to ablate the entire wafer equally.
- Spot mode : Multiple specified spots can be analyzed automatically. The spot information in KLARF format file can be imported.
- Straight mode : Straight line specified by two spots is analyzed with a specified width continuously.
- Block mode : It is similar to the straight mode, but a specified block size (e.g. 1 mm x 1 mm) is ablated for a specified time. After waiting for an interval time, the block next to the ablated block is ablated. This mode is better than the straight mode in terms of space resolution.
- Depth profile mode : The same spot with a specified size is ablated repeatedly to get depth profile information. An interval time can be set after ablating one layer, which gives better resolution of depth profile.
- Bevel mode : Multiple specified areas of wafer bevel can be analyzed automatically.

Expert\_LA model includes two load ports, which enables fully automated analysis of wafer samples by setting up a recipe for each wafer. A sample wafer is taken out from one of L/Ps, and it is recovered in another L/P after analysis of wafer.

LAGM software controls *MSAG\_DS* that has two special syringe pumps. One of syringe pumps contains a 1% HNO<sub>3</sub> solution and the other syringe pump contains a mixed standard solution in 1% HNO<sub>3</sub>, which are injected at 3 μL/min in total and the ratio of two syringe pumps are changed to make

## Specification and Environment

### Basic components

- Femto-second laser : 257 nm
- Maximum ablation frequency : 60 kHz
- Galvanometer mirror moving speed : 100 mm/sec
- X-Y-Z-θ stage
- Ejector
- GED-Q
- MSAG-DS
- Mass flow controllers
- LAGM intelligent software
- FFU (ISO class 3 environment)

### Option

- Expert\_LA that includes the followings:
  - FOUP loader : 2 set
  - Wafer transfer robot : 1 set

### Environment

Temperature : 15~30°C  
Humidity : 35~85%RH No condensation

Power supply, size and weight depends on the model. Please contact us.

IAS Inc.



2-2-1 Hinohonmachi, Hino, Tokyo 191-0011 Japan  
TEL: 042-589-5525 FAX: 042-589-5526  
E-Mail: [iasjapan@iasinc.jp](mailto:iasjapan@iasinc.jp) URL: <https://iasinc.jp>